

# Mepilex® Transfer

SafetaC®  
TECHNOLOGY



## The effective exudate transfer dressing with SafetaC® technology

- Minimizes pain and trauma at dressing changes<sup>1,2</sup>
- Transfers exudate away from the wound thereby minimizing the risk of maceration<sup>3</sup>
- Very convenient and conformable for difficult-to-dress wounds

### SafetaC® technology layer

- Reduces pain and trauma during dressing changes<sup>1</sup>
- Does not adhere to the moist wound bed, adheres only to dry skin
- Seals the wound margins and reduces risk of maceration<sup>3</sup>

### Thin polyurethane foam

- Transfers the exudate away from the wound into a secondary absorbent dressing
- Conforms well on difficult-to-dress wounds
- Ideal for large, exuding wounds



### SafetaC® technology. Less pain and less trauma.

SafetaC® is a patented adhesive technology that minimizes pain to patients and trauma to wounds. SafetaC® technology is available exclusively on Mölnlycke Health Care dressings, including Mepilex® wound care dressings, Mepitel®, Mepiform® and Mepitac®.

For more information visit [www.molnlyckeUSA.com](http://www.molnlyckeUSA.com)




Painful skin stripping with traditional adhesive




No skin stripping with SafetaC® technology


## How to use Mepilex® Transfer



Clean the wound area. Cut for customization, if desired. Remove the larger release film.



Mepilex® Transfer should overlap the wound bed by 1-2 cm in small product sizes and 5 cm in larger sizes onto the surrounding skin. Apply the adherent side to the wound. Remove the shorter release film and mold the dressing in place. Do not stretch.



Apply an appropriate absorbent secondary dressing over Mepilex® Transfer and fixate.

## How Mepilex® Transfer works

Mepilex® Transfer is a thin and highly conformable foam dressing that absorbs exudate and maintains a moist wound environment. The Safetac® layer seals around the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimizing the risk of maceration. As the Safetac® layer seals around the wound, the foam structure of Mepilex® Transfer allows the exudate to move vertically into a secondary absorbent dressing, maintaining a moist wound environment. The Safetac® layer ensures atraumatic dressing changes. Mepilex® Transfer can be cut to suit various wound shapes and locations.

## Benefits of Mepilex® Transfer

- Transfers exudate away from the wound
- Atraumatic to the wound and surrounding skin upon removal
- Conforms well to body contours
- May be cut to shape
- Minimizes the risk of maceration
- Promotes patient comfort during wear time and minimizes pain upon removal
- Suitable for large wounds and fragile tissue

## Note

If clinical signs of infection are present, the use of Mepilex® Transfer may be continued if proper infection treatment is initiated.

## Areas of use

Mepilex® Transfer is designed for a wide range of exuding wounds such as fungating tumours<sup>3</sup>, radiation-induced dermatitis, Epidermolysis Bullosa (EB), donor sites and venous leg ulcers. Mepilex® Transfer can also be used as a protective layer on non-exuding wounds and large areas with fragile skin.

**Wear time:** Up to 7 days



Radiation-induced dermatitis



Epidermolysis Bullosa (EB)



Fungating malignant wound

## Mepilex® Transfer Ordering Information†

Product Code	Size	Pcs/box	Pcs/case	HCPCS
294899	6" x 8" (15 x 20 cm)	5	40	A6210
294599	8" x 20" (20 x 50 cm)	4	24	A6211

† Packaged sterile in single packs

## References:

1. White R. A Multinational survey of the assessment of pain when removing dressings. Wounds UK 2008; Vol 4, No 1.
2. Dykes PJ et al. Effects of adhesive dressings on the stratum corneum of the skin. J Wound Care 2001; 10(2):7-10.
3. Dykes PJ. The effect of adhesive dressing edges on cutaneous irritancy and skin barrier function. J Wound Care 2007, 16(3):97-100.



Scan for more information

Socialize with us:



**We're here to help you, when you need us. Call your Mölnlycke® Health Care Representative or Regional Clinical Specialist.**

1-800-843-8497 | [www.molnlyckeUSA.com](http://www.molnlyckeUSA.com) | [info.woundcare@molnlycke.com](mailto:info.woundcare@molnlycke.com) | 5550 Peachtree Parkway, Suite 500, Norcross, GA 30092

The Mölnlycke Health Care, Mepilex, Mepitel, Mepiform, Mepitac, and Safetac, names and respective logos are registered trademarks of Mölnlycke Health Care AB. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2013. Mölnlycke Health Care AB. 1-800-843-8497.

MHC-2013-911 - 270012 05.13