



SIGVARIS

COMPREFLEX™

INELASTIC COMPRESSION WRAPS

**SIMPLE SOLUTIONS FOR
LYMPHOEDEMA
VENOUS INSUFFICIENCY
WOUND CARE**

COMPREFLEX BELOW KNEE



FEATURES

<21–46mmHg of gradient, short-stretch compression

Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture

Strap extenders may be applied if the patient exceeds standard size ranges (see page 14)

Combine with COMREFLEX THIGH COMPONENT for full leg compression

Available in black or beige

Easy care - machine wash and dry

Latex-free

INDICATIONS

Stage I or II Lymphoedema

Chronic Venous Insufficiency

Post-surgical Edema

Lipodermatosclerosis

INCLUDES

1 COMPREFLEX BK

1 COMPREBOOT

Set of ACCUTAB devices

Pair of BASIC liners

5 SIZES AVAILABLE

2 LENGTHS AVAILABLE

2 COLOURS AVAILABLE

BASIC LINER

Non-compressive




- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties



INCLUDES COMPREBOOT

- Ideal for mild-moderate edema/lymphoedema
- Features neoprene material with hook and loop tabs
- 4-way stretch for firm compression around ankle
- Provides moderate, therapeutic compression



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices

P.O. #	Account #	Contact	Date	
Phone	Patient Name	Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size	Length	Item #	Color
<input type="checkbox"/> <input type="checkbox"/>				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Measuring instructions on page 17

COMPREFLEX - BK SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

BLACK

REGULAR	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
TALL	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT


BEIGE

REGULAR	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
TALL	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

NOTE: If foot size doesn't match COMPREFLEX size, please order the COMPREFLEX NF (pg 4) with separate foot piece (pg 12).



COMPREFLEX NO FOOT



FEATURES

- <21–46mmHg of gradient, short-stretch compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 14)
- Good option for patients in need of a different foot option (size or boot)
- Easily worn under most standard footwear Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphoedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX BK
- Set of ACCUTAB devices
- Pair of BASIC liners

5 SIZES AVAILABLE
2 LENGTHS AVAILABLE
2 COLOURS AVAILABLE

BASIC LINER Non-compressive




- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN

Garment features a curved design to provide a conforming fit

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



EASY TO USE

Straps rolls back and attach to garment to prevent tangling

P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight

Company Name

Dropship Name & Address

PRODUCT INFORMATION


Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

NOTE: If foot garment is desired, please see page 12 for options and order form.

Measuring instructions on page 17

COMPREFLEX - NF SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm
BLACK					
REGULAR	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
TALL	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
BEIGE					
REGULAR	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
TALL	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT



COMPREFLEX LITE TRANSITION



FEATURES

- <21–46mmHg of gradient, short-stretch compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 14)
- Easily worn under most standard footwear Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphoedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX BK
- Set of ACCUTAB devices
- Pair of TRANSITION liners

5 SIZES AVAILABLE
2 LENGTHS AVAILABLE
2 COLOURS AVAILABLE

TRANSITION LINER

10–15mmHg in Foot and Ankle



- Soft, flat seam knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN

Garment features a curved design to provide a conforming fit

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



EASY TO USE

Straps roll back and attach to garment to prevent tangling

P.O. #	Account #	Contact	Date	
Phone	Patient Name	Age	Height	Weight

Company Name

Dropship Name & Address


PRODUCT INFORMATION

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Measuring instructions on page 17

COMPREFLEX LITE TRANSITION SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
A	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm
BLACK					
REGULAR	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
TALL	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
BEIGE					
REGULAR	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
TALL	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT



TRANSITION LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)
SMALL	19-26cm	29-39cm	<7½	<6½	1601-CK
MEDIUM	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK
LARGE	26-36cm	40-50cm	10-12	9-11	1603-CK
X-LARGE	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK
XX-LARGE	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK

COMPREFLEX LITE COMPLETE



FEATURES

- <21–46mmHg of gradient, short-stretch compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 14)
- Easily worn under most standard footwear Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphoedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX BK
- Set of ACCUTAB devices
- Pair of COMPLETE liners

5 SIZES AVAILABLE
2 LENGTHS AVAILABLE
2 COLOURS AVAILABLE

COMPLETE LINER

10–15mmHg in Foot, Ankle and Calf



- Comfortable knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

CONTOURED DESIGN


Garment features a curved design to provide a conforming fit

EASY TO USE

Straps rolls back and attach to garment to prevent tangling

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices



P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION


Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Leg	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Measuring instructions on page 17

COMPREFLEX LITE COMPLETE SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
B	24–34cm	29–39cm	34–44cm	39–49cm	44–55cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
BLACK					
REGULAR	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
TALL	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT
BEIGE					
REGULAR	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
TALL	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT



COMPLETE LINER SIZING CHART




SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)
SMALL	19–26cm	29–39cm	<7½	<6½	1601-CK-F1
MEDIUM	21–30cm	34–44cm	7.5–9.5	6.5–8.5	1602-CK-F1
LARGE	26–36cm	40–50cm	10–12	9–11	1603-CK-F1
X-LARGE	31–41cm	46–56cm	12.5–15	11.5–14	1604-CK-F1
XX-LARGE	36–46cm	54–64cm	12.5–15	11.5–14	1605-CK-F1

EASILY AND ACCURATELY SET COMPRESSION WITH ACCUTAB




EASY TO USE

1. Gently lay the COMPREFLEX straps in place without stretching the fabric
2. Attach an ACCUTAB at the end of each strap
3. Pull each strap to the prescribed dosage indicated by the lines on the ACCUTAB:
<21, 23–32, or 34–46mmHg



PROVEN ACCURACY FOR RESTING PRESSURE

ACCUTAB RESULTS: 20–30mmHg



* ACCUTAB has been rigorously tested according to established guidelines in the International Compression Club Consensus documents 1, 2. Testing was conducted using PicoPress®, which is recognized for its high accuracy and reproducibility 3.

EFFECTIVE WORKING PRESSURE


- Inelastic compression reduces leg edema and improves venous blood flow 4, 5
- Inelastic compression can increase arterial perfusion in patients with mild to moderate Peripheral Arterial Disease 4, 5

WORKING PRESSURE

Walking, Stair Climbing, Sit to Stand, Toe Raises

RESTING PRESSURE

Donning
COMPREFLEX



RESTING PRESSURE

Laying Down
Sitting Still

SAFE, PRECISE RESTING PRESSURE

- Inelastic compression can be used safely for patients with mild to moderate Peripheral Arterial Disease 4, 5 *
- ACCUTAB delivers precision dosages for controlled initial resting pressure when donning SIGVARIS COMPREFLEX garment

COMPREFLEX LEG WRAP

HIGH QUALITY

Each garment is double stitched to ensure durability and quality standards

CONTOURED DESIGN

Garment features a curved design to provide a conforming fit

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance

EASY TO USE


Straps roll back and attach to garment to prevent tangling

COMFORTABLE MATERIALS

Soft, conforming Breathe-O-Prene® provides a comfortable fit and wicks away moisture

TWO Colour OPTIONS

The COMPREFLEX is available in black or beige



ACCUTAB

Each COMPREFLEX lower leg wrap includes a set of ACCUTAB compression measuring devices

COMPREFLEX THIGH COMPONENT



FEATURES

- Includes ACCUTAB for easy, accurate compression
- Soft, conforming Breathe-O-Prene provides a comfortable fit and wicks away moisture
- Strap extenders may be applied if the patient exceeds standard size ranges (see page 14)
- Combine with COMPREFLEX BK for full leg compression
- 2 KNEE options available
- Available in black or beige
- Easy care - machine wash and dry
- Latex-free

INDICATIONS

- Stage I or II Lymphoedema
- Chronic Venous Insufficiency
- Post-surgical Edema

INCLUDES

- 1 COMPREFLEX TC
- Set of ACCUTAB devices
- Belt & Hip Attachment
- Pair of non-compressive cotton liners

4 SIZES AVAILABLE
2 LENGTHS AVAILABLE **L + R**
2 COLOURS AVAILABLE



BELT & HIP ATTACHMENT

These items are included to ensure a comfortable fit, and help secure garment in place

EASY TO USE

Straps rolls back and attach to garment to prevent tangling

STRETCH PANEL

The front stretch panel secures garment in place for easy application and high compliance

KNEE OPTIONS *(additional charge)*



COMPREFLEX KNEE

Closed Patella

Soft, conforming compression with straps that wrap around the back of the knee



COMPREFLEX KNEE

Open Patella

Soft, conforming compression with straps that wrap around the back of the knee

3 SIZES AVAILABLE
2 COLOURS AVAILABLE

P.O. #	Account #	Contact	Date
Phone	Patient Name	Age	Height

Company Name

Dropship Name & Address

PRODUCT INFORMATION

Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Leg	Size/Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

KNEE STYLE	Item #	Size	Color
<input type="checkbox"/> OPEN PATELLA <input type="checkbox"/> CLOSED PATELLA			<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

Measuring instructions on page 17

COMPREFLEX - TC SIZING CHART (D-F)

	SMALL	MEDIUM	LARGE	XLARGE
F	48-58cm	56-66cm	64-74cm	74-84cm
E	43-53cm	51-61cm	58-68cm	68-78cm
D	38-48cm	46-56cm	53-63cm	63-73cm


LEFT					
BLACK	REGULAR	1401-TCR-L	1402-TCR-L	1403-TCR-L	1404-TCR-L
	TALL	1401-TCT-L	1402-TCT-L	1403-TCT-L	1404-TCT-L
RIGHT					
BLACK	REGULAR	1401-TCR-R	1402-TCR-R	1403-TCR-R	1404-TCR-R
	TALL	1401-TCT-R	1402-TCT-R	1403-TCT-R	1404-TCT-R
LEFT					
BEIGE	REGULAR	1411-TCR-L	1412-TCR-L	1413-TCR-L	1414-TCR-L
	TALL	1411-TCT-L	1412-TCT-L	1413-TCT-L	1414-TCT-L
RIGHT					
BEIGE	REGULAR	1411-TCR-R	1412-TCR-R	1413-TCR-R	1414-TCR-R
	TALL	1411-TCT-R	1412-TCT-R	1413-TCT-R	1414-TCT-R

COMPREKNEE SIZING CHART

	SMALL	MEDIUM	LRG / X-LRG
D	38-48	46-56	53-73
C1	29-39	34-44	39-65

CLOSED PATELLA			
BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP

OPEN PATELLA			
BLACK	1101-KP-OP	1102-KP-OP	1103-KP-OP
BEIGE	1111-KP-OP	1112-KP-OP	1113-KP-OP



FOOT COMPRESSION

COMPRESSBOOT LITE



2 SIZES AVAILABLE
2 LENGTHS AVAILABLE
COMPRESSBOOT

10–15mmHg

- Ideal for mild edema/lymphoedema
- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design (may be trimmed to fit)
- Breathable stretch fabric conforms to ankle
- Provides mild, therapeutic compression

COMPRESSBOOT LITE PRODUCT CODES

	SMALL/ MEDIUM / LARGE	XLRG / XX LRG
REGULAR	1803-BLR	1805-BLR
LONG	1803-BLL	1805-BLL

30–40mmHg

- Ideal for mild-moderate edema/lymphoedema
- 4-way stretch for firm compression around ankle
- Features neoprene material
- Provides moderate, therapeutic compression

with hook and loop tabs

COMPRESSBOOT PRODUCT CODES

	SMALL	MEDIUM/LARGE	XLRG / XX LRG
BLACK			
REGULAR	1801-BTR	1803-BTR	1805-BTR
LONG	1801-BTL	1803-BTL	1805-BTL
BEIGE			
REGULAR	1811-BTR	1813-BTR	1815-BTR
LONG	1811-BTL	1813-BTL	1815-BTL



3 SIZES AVAILABLE
2 LENGTHS AVAILABLE
2 COLOURS AVAILABLE
COMPRESSBOOT PLUS

30–40mmHg

- Ideal for chronic/acute edema/lymphoedema
- Chipped foam lining
- Removable gripper patches provide traction
- Straps are removable for custom positioning or replacement
- Provides moderate, therapeutic compression

COMPRESSBOOT PLUS PRODUCT CODES

	SMALL	MEDIUM / LARGE	XLRG / XX LRG
REGULAR	1801-BPR	1803-BPR	1805-BPR
LONG	1801-BPL	1803-BPL	1805-BPL



3 SIZES AVAILABLE
2 LENGTHS AVAILABLE


P.O. #	Account #	Contact	Date		
Phone	Patient Name		Age	Height	Weight
Company Name					
Dropship Name & Address					

PRODUCT INFORMATION	
Foot <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	
Size	<input type="checkbox"/> COMPREBOOT LITE
Length	<input type="checkbox"/> COMPREBOOT
Item #	<input type="checkbox"/> COMPREBOOT PLUS

MEASURING INSTRUCTIONS

- I:** Measure circumference of ankle bend and heel.
- J:** Measure circumference across metatarsal heads.
- K:** Measure length on the medial side of the foot from the heel to the 1st metatarsal head.

CIRCUMFERENCE	
LEFT	RIGHT
_____	_____ I
_____	_____ J
_____	_____ K



FOOT SIZING CHART (I-K)						
	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

ACCESSORIES & LINERS



COMPREFLEX STRAP EXTENDER

Adds 10 cm of adjustable range to a single strap.

1 SET

BLACK	1401-SE1
BEIGE	1411-SE1



COMPREFLEX KNEE

Compression garment with either open or closed patella.

	SMALL	MEDIUM	LRG / X-LRG
D	38-48	46-56	53-73
C1	29-39	34-44	39-65

CLOSED PATELLA

BLACK	1101-KP	1102-KP	1103-KP
BEIGE	1111-KP	1112-KP	1113-KP

OPEN PATELLA

BLACK	1101-KP-OP	1102-KP-OP	1103-KP-OP
BEIGE	1111-KP-OP	1112-KP-OP	1113-KP-OP



BASIC LINER


Non-compressive

- Easy to apply knee-high liners
- *drirelease* moisture wicking and odor control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

Size: One size fits all

PAIR

1700-NK1



ACCUTAB

Compression measuring device design to use with COMPREFLEX leg wraps.

SIZE 1 SET OF 5

SMALL	1491-TB1
MEDIUM	1492-TB1
LARGE	1493-TB1
X-LARGE	1494-TB1
XX-LARGE	1495-TB1



COMPLETE LINER

10-15mmHg in Foot, Ankle & Calf

- Comfortable knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties



TRANSITION LINER

10-15mmHg in Foot & Ankle

- Soft, flat seam knee-high liners
- *drirelease* moisture wicking and odor control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

COMPLETE LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)
SMALL	19-26cm	29-39cm	<7½	<6½	1601-CK-F1
MEDIUM	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK-F1
LARGE	26-36cm	40-50cm	10-12	9-11	1603-CK-F1
X-LARGE	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK-F1
XX-LARGE	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK-F1

TRANSITION LINER SIZING CHART

SIZE	ANKLE CIRC. (A)	CALF CIRC. (C)	WOMEN'S SHOE	MEN'S SHOE	PAIR (BLACK)
SMALL	19-26cm	29-39cm	<7½	<6½	1601-CK
MEDIUM	21-30cm	34-44cm	7.5-9.5	6.5-8.5	1602-CK
LARGE	26-36cm	40-50cm	10-12	9-11	1603-CK
X-LARGE	31-41cm	46-56cm	12.5-15	11.5-14	1604-CK
XX-LARGE	36-46cm	54-64cm	12.5-15	11.5-14	1605-CK

Follow contour of the arm/leg by placing the measuring tape against the skin. Measure in centimeters.

BELOW KNEE

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 15cm, 25cm and 30cm dots. Record in the corresponding lines A, B, C, C1(Left/Right).

LENGTH

Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**, and record in box G.

FOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

THIGH COMPONENT

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the lateral aspect **Top of Patella** and continue to the 15cm, and 23cm dots. Record in the corresponding lines D, E, F (Left/Right).

LENGTH (COMPREFLEX only)

Measure length from the lateral aspect **Top of Patella** to lateral aspect **Gluteal Fold**, and record in box G1.

THIGH HIGH

CIRCUMFERENCES

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone) and continue to the 15cm, 25cm and 30cm dots. Record in the corresponding lines A, B, C, C1(Left/Right). Measure circumferences at the 5cm dot from the lateral aspect **Top of Patella** and continue to the 15cm, and 23cm dots. Record in the corresponding lines D, E, F (Left/Right).

LENGTH

Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Gluteal Fold**, and record in box H.

FOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

*For further information and fitting contact our SIGVARIS Specialist:
0800 SIGVARIS (744 827) | sigvaris@uslmedical.co.nz | www.uslmedical.co.nz*


SIGVARIS COMPRESSION STOCKING & SLEEVES

TRADITIONAL - Class 2 (503) (23-32 mmHg) Class 3 (504) (34-46 mmHg) Class 4 (505) minimum mmHg)

When a high degree of safety and effectiveness is required.

- Maximum pressure stability throughout the day
- Excellent strength-stretch ratio
- Great medical benefits for safe and successful treatment
- High wearing comfort thanks to special elasticity and good air permeability
- With refined natural rubber
- Wide range of sizes available

When it comes down to high medical efficacy, TRADITIONAL is the reliable classic. Day after day, this product gives your legs a secure feeling, ensuring top-quality medical compression. As TRADITIONAL stockings are particularly elastic and breathable, they are extremely comfortable to wear. The neutral skin colour matches all clothing and combines perfectly with delicate stockings.



COTTON - Class 1 (18-21 mmHg) Class 2 (23-32 mmHg) Class 3 (34-46 mmHg)

COTTON stockings make your legs feel light and provide exceptional comfort. The fine, processed cotton yarn is lovely and soft, and is in direct contact with the skin. Their excellent thermo-regulation makes COTTON stockings extremely comfortable to wear in both cold and warm weather.

- Optimal fit and maximum comfort
- Exclusive cotton yarn inside the stocking (SWISS COTTON Premium cotton fibres)
- Ideal for men and women with sensitive skin
- Excellent thermo-regulation
- Easy to slip on and off
- SWISS MADE quality seal
- Wide range of sizes available
- Flesh and Black colours available




MAGIC - Class 1 (18-21 mmHg) Class 2 (23-32 mmHg)

The women's compression stocking for every requirement.

- A delicate, slightly transparent stocking with a natural, elegant look
- Perfect fit and exceptional degree of comfort
- Simple handling
- Provides optimum leg support
- Wide range of sizes available
- Flesh and Black colours available

MAGIC, the versatile compression stocking, is exceptionally easy to handle and has a perfect fit. This has been confirmed by intensive tests carried out with compression stocking wearers under scientific conditions. Their natural look and subtle transparency also appeal to the senses.



ARM SLEEVES - Traditional and Advance

Advance 1 + 11

Silver ion enriched polymide microfibre has antibacterial properties and prevents odours and skin irritations.

- Bacteriostatic polymide microfibre
- Thermoregulating fibres
- Sensinova griptop or soft top
- With or without mitten



Traditional

For maximum pressure stability throughout the day.

- Griptop or arm strap
- With or without mitten
- Mitten available separately



SIGVARIS

COMPREFLEX™

INELASTIC COMPRESSION WRAPS



Phone

0800 (SIGVARIS) 744 827

Email

sigvaris@uslmedical.co.nz

* Safe compression dosage should be prescribed by a physician based on individual patient medical diagnoses and needs. Compression should not be applied on patients with severe Peripheral Arterial Disease as tissue necrosis, limb loss, or death may result. When using ACCUTAB, do not apply more mmHg of resting compression pressure than prescribed by your physician. Immediately remove your compression garment and consult with your physician if you experience leg pain, numbness or tingling, coolness, or any other adverse symptoms.

(1) Partsch, Hugo, et al. "Measurement of lower leg compression in vivo: recommendations for the performance of measurements of interface pressure and stiffness." *Dermatologic Surgery* 32.2 (2006): 224–233.

(2) <http://www.tagungsmanagement.org/comp/index.php/icc-consensus>

(3) Mosti, Giovanni, and Susanna Rossari. "The importance of measuring sub bandage pressure and presentation of new measuring device." *Acta Vulnol* 6 (2008): 31–36.

(4) Mosti, Giovanni, Maria Letizia Iabichella, and Hugo Partsch. "Compression therapy in mixed ulcers increases venous output and arterial perfusion." *Journal of vascular surgery* 55.1 (2012): 122–128.

(5) Mosti, Giovanni. "Compression in mixed ulcers: venous side." *Phlebology* 29.1_suppl (2014): 13–17.

SIGVARIS, COMPREFLEX, Breathe-O-Prene, and ACCUTAB are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2017 SIGVARIS, Inc. 2017 *driverelease*. *driverelease* is an Optimizer Brands company. All Rights Reserved.

LIT-COMPREFLEXBOOK 09.17

Distributed by



Order by freephone:
Auckland: 0800 658 814
Christchurch: 0800 804 546

Order by freefax:
Auckland: 0800 830 660
Christchurch: 0800 804 547

Auckland Head Office
494 Rosebank Road, Avondale, Auckland 1026
Email: customerservice@uslmedical.co.nz
www.uslmedical.co.nz