

COMPREFLEX REDUCE LOWER LEG

P.O. #	Account #	Contact	Date
Phone	Patient Name	Age	Height Weight
Company Name			
Dropship Name & Address			

PRODUCT INFORMATION

COMPREFLEX REDUCE BK	Size	Length	Item #	Color
☐ LEFT ☐ RIGHT				☐ BLACK ☐ BEIGE
COMPREFLEX REDUCE BOOT	Size		Item #	Color
☐ LEFT ☐ RIGHT				☐ BLACK ☐ BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS

COMPREFLEX REDUCE BELOW KNEE

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 25cm dot. Record in the corresponding lines A and C (Left/Right). Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**, and record in box G.

COMPREFLEX REDUCE BOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

COMPREFLEX REDUCE BK SIZING CHART

	RANGE 1	RANGE 3
C	<65cm	<90cm
A	<55cm	<75cm
BLACK		
REGULAR	2401-BKR	2403-BKR
TALL	2401-BKT	2403-BKT
BEIGE		
REGULAR	2411-BKR	2413-BKR
TALL	2411-BKT	2413-BKT

COMPREFLEX REDUCE BOOT SIZING CHART

	RANGE 1	RANGE 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
BLACK	2801-BTL	2803-BTL
BEIGE	2811-BTL	2813-BTL

CIRCUMFERENCE

LEFT RIGHT

