

## Pain Relief

**Vectored compression:** Delivers targeted compression to the pain site

**Integrated hot/cold gel pad:** Easily slides into the lumbar pouch **Adjustable lumbar panel:** Shifts right or left for precise focus of the compression

Clinical validation: A recent study using the bilateral cinching mechanism of the Vector showed "Patient's average VAS pain score dropped by 41% immediately upon donning the device ... "

\*"Prospective Study of Lumbar Back Brace In an Interventional Pain Practice." Rosenthal, Richard M.D. & Spencer, Shawn, BS. Practical Pain Management, (July/August 2011)

## Comfort

**Breathable spacer fabric:** Allows heat dissipation and quick evaporation for prolonged use

**Contoured waist belt:** Fits comfortably around the rib cage and iliac crest - naturally conforms to all body shapes

Lightweight: Weighs less than 2 pounds

## Ease of Use

**Bilateral cinching:** Bilateral cinching mechanism doubles the compression provided with equal pull force

**Mechanically-advantaged pulleys:** Pulleys reduce the force needed to apply compression

**Ergonomic hand pulls:** Substantial hand pulls are easy to grasp, even for arthritic patients


Easy-release Velcro: Provides firm grip and easy release

## Indications

Back Sprains & Strains Herniated Discs

- · Degenerative Disc Disease
- Lumbago





Bilateral Cinching



Mechanically-advantaged pulleys



Ergonomic hand pulls



Hot/Cold Gel Pad



Contoured Fit



Breathable Fabric

