ProShield®

Mask Range




THERAPIES. HAND IN HAND.



PROSHIELD® N-95

- Meets NIOSH criteria for protection against TB
- Secure full facial fit
- Easy donning
- Conformable nose piece
- Durable and comfortable headbands
- Optimum protection against blood and body fluid splash (160mmHg)







With a **Bacterial Filtration Efficiency (BFE)** of greater than 99% for media greater than 3 microns this mask provides superior protection from potential infection.¹



The super high **Particulate Filtration Efficiency (PFE)**material filters >99% of particles
greater than 0.1 microns. When
securely fitted, the mask helps to
provide an effective barrier to
particulate matter from laser plumes
and electrocautery procedures.²



The N95 mask is NIOSH approved as a N95 particulate filter respirator. It meets or exceeds the standard performance criteria demanded by NIOSH for the management of Tuberculosis (TB).³



The Fluid Resistant qualities of the mask provide protection against fluid strikethrough.⁴

File Notes:

- BFE (Bacterial Filtration Efficiency %) is an independent test conducted by Nelson Laboratories that measures the effectiveness of face mask materials to filter out bacteria.
- 2. PFE (Particulate Filtration Efficiency) measures the effectiveness of the mask material to filter particulate matter.
- 3. NIOSH National Institute of Occupational Safety and Health (USA).
- 4. Fluid resistance is a standard test for resistance of medical face masks to penetration by synthetic blood.

ORDERING INFORMATION

Ref-no.	Description	Items per unit	Shipper Qty
72509-09	ProShield® N-95 Small	bx/50	6 bxs
72509-10	ProShield® N-95 Medium	bx/50	6 bxs

BSN medical Australia & New Zealand

L6, 211 Wellington Road, Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au www.bsnmedical.co.nz Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand T 0508 276 111 F 0508 998 830

