## RHS 402 // IPA SURFACE DISINFECTION WIPES




Reynard's IPA wipes contain 70% Isopropyl Alcohol, proven to be effective against 99.9% of common germs. Presented in the familiar canister format they are highly portable, convenient and practical.

A flip top lid prevents the wipes from drying out. They have been tested and proven to be effective against a wide range of bacteria: with a kill rate of 99.9% of common germs.

Directions: To ensure wipes are saturated with the fluid, tip the canister end to end. Remove the wipes from the canister as needed and disinfect surfaces. Dispose of the used wipes in a suitable container.





| Features  | Benefits   |
|---|--|
| <ul> <li>Large quick drying wipes</li> <li>70% Isopropyl Alcohol</li> <li>TGA Listed ARTG 164342</li> <li>NZFSA Approved C44 (non dairy)</li> </ul> | <ul> <li>Complete disinfection</li> <li>Effective against MRSA &amp; E.coli</li> </ul> |

## **Applications**

- ✓ Hospitals
- ✓ Rest homes
- ✓ Factories